GUIDING PHARMACEUTICAL MANUFACTURERS IN CLEANROOM GLOVE SELECTION: BALANCING PERFORMANCE, COMPLIANCE, AND SUSTAINABILITY.

Walid El Azab

QP Pro Services & QPM Consulting
Co-founder & Senior Consultant

Trusted **for** **Continuity.**

Guiding Pharmaceutical Manufacturers in Cleanroom Glove Selection: Balancing Performance, Compliance, and Sustainability.

Walid El Azab – QP Pro Services & QPM Consulting
Co-founder & Senior Consultant

Introduction

In pharmaceutical manufacturing, cleanroom gloves are indispensable for maintaining contamination control, ensuring product integrity, and safeguarding operator health. Whether in sterile or non-sterile or laboratory environments, the selection of cleanroom gloves is a critical decision that hinges on various factors, including glove safety, functionality, and performance for the intended use. Another factor weighing the balance is the growing emphasis on sustainability, including Environmental, Social, and Governance (ESG) factors.

This article provides a comprehensive framework to guide pharmaceutical manufacturers in making informed decisions about cleanroom glove selection. It emphasizes the importance of balancing performance, compliance, and sustainability. By integrating regulatory requirements and ESG priorities, the article aims to empower stakeholders to align glove selection with corporate sustainability objectives while supporting operational compliance. Finally, it presents a decision tree as a practical tool for glove selection.

Key Considerations in Cleanroom Glove Selection

Cleanroom gloves used in pharmaceutical applications must comply with strict regulatory and quality standards to ensure safety, functionality, and performance. In Europe, these gloves fall under EU Regulation 2016/425, which classifies Personal Protective Equipment (PPE) into three categories (1):

- Category I: PPE for minimal risks.
- Category II: PPE for moderate risks.
- Category III: PPE designed to protect against serious or irreversible health hazards.

Given the potential exposure to chemical and biological hazards in pharmaceutical environments, cleanroom and laboratory gloves are expected to be certified as Category III PPE to ensure adequate protection.

Glove manufacturers must meet essential health and safety requirements, undergo conformity assessments, and affix CE markings to demonstrate compliance. Responsibility for product safety spans the supply chain, including manufacturers, representatives, importers, and distributors.

Cleanroom gloves must comply with EU and international standards for safety and performance in pharmaceutical environments. General requirements like fit, comfort, and usability are evaluated under ISO 21420:2020. Durability is crucial to ensure prolonged use and protective integrity; tests include:

- Pinhole defects: ASTM D5151-19
- Tensile strength and elongation: ASTM D412-16
- Cleanliness: IEST-RP-CC005.4, ASTM D6124-06
- Chemical resistance and barrier efficacy: ISO 374 series (ISO 374-1:2024, ISO 374-2:2019, ISO 374-4:2019, EN 16523-1:2015+A1:2018)
- Microbial barrier integrity: ASTM F1671/F1671M-22, ISO 374-5:2024

Compliance with additional standards or test methods may be required based on the end user's needs. For example, gloves intended for aseptic areas may need to meet ISO 11137 series for sterilization, while gloves used in chemotherapy drug manufacturing may need to comply with USP <800> or be tested according to ASTM D6978-05:2023 for resistance to cytotoxic drugs.

Factors to Consider for Glove Selection

The industry has transitioned from PVC and latex gloves to nitrile gloves due to their optimal hypoallergenic properties, chemical resistance, flexibility, and cost (2-5). Latex gloves, though highly elastic, carry risks of Type I hypersensitivity and batch-to-batch variability due to their natural origin (6-8). Vinyl gloves, on the other hand, are less durable and lack the chemical resistance (2-5). Nitrile gloves offer a balanced combination of durability, flexibility, safety, and cost, making them the preferred choice (2-8).

The glove's safety, functionality, and performance depend on the manufacturing process, which ensures quality, cleanliness, and compliance with relevant standards. Understanding these manufacturing steps is crucial for informed decision-making:

1. **Glove Forming Preparation:** The process begins with cleaning hand molds by immersing them sequentially in acid, water, and alkali solutions. This step removes residual impurities to ensure the purity of the finished gloves.
2. **Immersion Forming and Vulcanization:** Cleaned molds are coated with a coagulant mixture and immersed in latex or nitrile formulations. Chemical accelerators like thiurams or dithiocarbamates facilitate vulcanization by cross-linking polymer chains. Vulcanization enhances elasticity, durability, and tear resistance. However, these accelerators may cause Type IV hypersensitivity, driving the industry's shift toward accelerator-free nitrile gloves (8).
3. **Chlorination and Surface Treatment:** After vulcanization, gloves undergo chlorination to reduce tackiness and improve wearability and comfort. Secondary chlorination, often performed offline, ensures optimal cleanliness and smoothness by targeting residual contaminants on the external surface.
4. **Cleanroom Processing:** Rigorous washing cycles (e.g. with filtered deionized water) are applied, followed by drying with HEPA-filtered air. Gloves are packaged in cleanroom environments, and gloves labeled as sterile undergo sterilization using gamma irradiation or ethylene oxide.

A. Glove safety factors:

1. **Compatibility with processes:** Glove materials must resist chemical degradation, frequent disinfection, and mechanical stress to ensure prolonged usability and reliability. Selecting the right glove material requires alignment with specific manufacturing demands, including chemical exposure, mechanical stress, and contamination control requirements.
 - **Process-Specific Needs:** For instance, gloves used in chemotherapy drug manufacturing must comply with ASTM D6978-05:2023 and USP <800> to ensure adequate chemical barrier properties. Gloves used in wet operations must resist permeation by chemical cleaning agents or solutions. Conversely, gloves for dry operations should prioritize electrostatic discharge properties to prevent static-related risks.
 - **Aseptic Cleanrooms:** Using sterile gloves is mandatory to maintain contamination control in aseptic cleanrooms. Gloves are often worn continuously for up to 4 hours and require frequent sanitization—typically every 5 to 10 minutes—due to contact with critical surfaces and materials. This frequent disinfection process, primarily involving 70% isopropyl alcohol (IPA), is essential for effective microbial control (10). However, repeated

exposure to IPA can degrade glove integrity materials over time, increasing the risk of micro-perforations and potentially compromising contamination control (11-13). To mitigate this risk, implementing regular glove replacement schedules is crucial for preserving glove integrity, maintaining performance, or ensuring aseptic conditions within cleanroom environments (14).

2. **Allergen Mitigation:** Latex-free gloves, especially accelerator-free nitrile gloves, effectively reduce the risks of Type I and Type IV hypersensitivity (6-9). Accelerator-free nitrile gloves help prevent allergic reactions, which is consistent with the EU REACH Regulation (EC 1907/2006) that restricts hazardous substances.
3. **Cleanliness for Cleanrooms:** Cleanroom gloves must adhere to the cleanroom cleanliness and intended use:
 - **Low Particle Shedding and Residue-Free Standards:** The gloves should meet the specific cleanliness based on the process and environment where used as such any gloves used in cleanroom should not generate particles or contaminant that could contaminate the cleanroom particle cleanliness where limits are described in the ISO 14644-1:2015. Therefore, glove manufacturer would provide evidence of particle counts or cleanliness proprieties using test method as per IEST-RP-CC005.4 for cleanroom gloves or as per ASTM D6124-06(2022) for laboratory gloves to ensure that residual powder levels remain under 2 mg, meeting industry benchmarks for residue-free performance.
 - **Microbial Barrier Properties:** Gloves used in operations involving contact with human organs, materials, other blood derivatives or other biological agents must demonstrate penetration resistance.

B. Glove comfort factors:

Operator comfort is a critical factor influencing performance and adherence to standard operating procedures in cleanroom environments.

1. **Ergonomic Design:** Gloves with a thin fingertip and palm design, along with high flexibility (tensile strength), enhance tactile sensitivity; a crucial feature for aseptic or complex/accurate operations requiring delicate material handling (3,4).
2. **Proper Fit:** Gloves must be appropriately sized to maintain dexterity; considerations include:
 - Texture: Textured fingertips and palms enhance grip, improving control over tools and materials.
 - Grip: Smooth gloves facilitate easy double-gloving, while tacky surfaces improve secure handling of delicate objects
3. **Double Gloving:** In aseptic environments, double gloving minimizes contamination risks (10). Standards protocol involves tucking sterile garment sleeves into a second pair of sterile gloves worn over the first pair donned during gowning (10). Smooth gloves with non-stick properties facilitate donning and doffing without compromising the cleanroom specifications.
4. **Regular Glove Changes:** Gloves should be replaced after critical interventions, loss of integrity, or extended use to maintain integrity or aseptic conditions. Ambidextrous glove designs simplify frequent changes, reducing disruptions in workflow and minimizing the risk of contamination.

C. Gloves Cost factors

The Total Cost of Ownership (TCO) for cleanroom gloves extends beyond initial purchase costs to highlight their long-term value. Pharmaceutical companies are implementing ESG principles into their operations to meet regulatory and social requirements. Amgen has constructed a sustainable manufacturing

facility in Ohio aligned with its 2027 environmental sustainability goals (15). Novartis is committed to achieving full carbon neutrality by 2030, meeting objectives outlined in the EU Green Deal and U.S. climate targets (16). Novo Nordisk has set goals to achieve zero environmental impact by 2045 (17). As such, (ESG) considerations play a vital role in cleanroom glove selection:

1. **Environmental perspective**, Effective disposal methods, such as recycling or downcycling, align with the EU Waste Framework Directive (2008/98/EC) and the U.S. Resource Conservation and Recovery Act (RCRA), minimizing environmental impact (18-20). Eco-friendly production methods, such as energy-efficient processes and CO₂ emissions reduction and other methods, comply with the EU and U.S. climate initiatives (21-23). For example, glove manufacturers adopting renewable energy technologies and water-saving innovations contribute to regional and global sustainability goals.
2. **Social considerations** are equally important. Compliance with the EU's Corporate Sustainability Reporting Directive (CSRD) and the U.S. Fair Labor Standards Act (FLSA) ensures ethical sourcing and fair labor practices (22, 24). Initiatives like the Responsible Glove Alliance (RGA) further address labor rights within the glove manufacturing industry (25). Health and safety are also paramount; the adoption of accelerator-free gloves mitigates the risks of chemical exposure and allergic reactions, supporting workplace safety (26).
3. **Governance** is another critical dimension of ESG. Supplier transparency is essential, with the EU's Corporate Sustainability Due Diligence Directive (CSDDD) requiring companies to ensure their suppliers adhere to ethical and environmental standards (27). Similarly, the U.S. SEC Climate Disclosure Rules mandate transparency in ESG impacts, fostering accountability (28).

End-User Decision Tree for Glove Selections:

A decision tree designed as a table-based (Table 1) selection process integrates compliance, operational needs, material properties, and ESG factors, enabling systematic comparison of glove options across suppliers.

Table 1: Decision tree represented as a table-based selection process

Question	Criteria/Standard	Response	Action
Step 1: Regulatory and Compliance Evaluation			
Does the glove meet regulatory requirements?	CE marking, EU Regulation 2016/425, OSHA standards, ISO 374-1/5, ASTM D6978, etc.	Yes / No	Exclude gloves that do not meet standards.
Does the glove pass essential testing?	ASTM D412 (tensile strength), ASTM F1671 (microbial barrier), ISO 374 (chemical resistance), etc.	Yes / No	Exclude gloves failing any key tests.
Step 2: Pharmaceutical Operational Needs			
Is the glove sterile for aseptic applications?	ISO 11137 sterilization, double-gloving compatibility, resistance to IPA degradation.	Yes / No	Exclude non-compliant gloves.
Does the glove resistant to chemicals hazardous?	ASTM D6978 compliance for chemotherapy drugs or ISO 374-1 for chemical permeation resistance.	Yes / No	Reject gloves unsuitable for chemical hazards.
Does the glove resistant to biological hazardous?	Microbial barrier integrity maintained under ASTM F1671/F1671M-22 or ISO 374-5:2024	Yes / No	Reject gloves unsuitable for biological hazards.
Does the glove meet cleanroom cleanliness?	IEST-RP-CC005.4 (cleanliness proprieties shared) in line with ISO 14644 where the	Yes / No	Exclude gloves failing cleanliness

	glove will be used (particle shedding in cleanroom) or ASTM D6124 (residual powder for laboratory gloves in boxes). Including rigorous washing and packaging in clean environments.		tests or glove with cleanliness proprieties not in line with the cleanroom specification as per ISO 14644.
Step 3: Material and Design			
Is the material suitable for the operation?	Nitrile (preferred for durability and hypoallergenic properties), latex (check allergy risks), vinyl (low-risk areas).	Yes / No	Reject gloves made from unsuitable materials.
Is the glove ergonomic and user-friendly?	Thin fingertips, textured surfaces, flexibility for prolonged use, proper sizing for dexterity.	Yes / No	Exclude gloves lacking ergonomic design.
Is the glove accelerator-free?	Compliance with EU REACH Regulation to reduce Type IV hypersensitivity risk.	Yes / No	Reject gloves with harmful accelerators.
Step 4: ESG Considerations			
Is the glove sustainably manufactured?	Biodegradable materials, energy-efficient production, reduced CO ₂ emissions.	Yes / No	Exclude gloves not meeting ESG criteria.
Does the supplier follow ethical labor practices?	Compliance with CSRD, FLSA, and Responsible Glove Alliance standards.	Yes / No	Reject suppliers not adhering to ethics.
Is the supply chain transparent?	ISO 14001 (environmental management), SEC Climate Disclosure compliance, documented ESG metrics.	Yes / No	Reject suppliers lacking transparency.
Step 5: Cost and Performance Comparison			
Does the glove have a favorable TCO?	Durable materials, long lifecycle, minimal need for replacements.	Yes / No	Reject gloves with high replacement costs.
Does the glove perform comparably across suppliers?	Evaluate supplier KPIs: chemical resistance, durability, cleanliness, ESG adherence.	Yes / No	Select gloves offering the best balance.

To use this table, begin by populating the responses for each question. For every glove option, record the response as either "Yes" or "No." Next, filter out non-compliant gloves by excluding any that received a "No" response at critical decision points. Finally, compare the remaining gloves, assessing them based on their performance, cost, and alignment with ESG criteria to make the final selection.

Conclusion

Selecting cleanroom gloves for pharmaceutical manufacturing requires balancing safety, functionality, regulatory compliance, and sustainability. Nitrile gloves emerge as the preferred choice due to their hypoallergenic properties, chemical resistance, and durability while meeting stringent regulatory standards. Compliance with ISO and ASTM guidelines ensures mechanical strength, microbial barrier effectiveness, and cleanliness. The adoption of accelerator-free nitrile gloves further minimizes allergenic risks, aligning with EU REACH regulations.

By integrating ESG principles, such as sustainable manufacturing practices and ethical sourcing, manufacturers can enhance both regulatory compliance and environmental stewardship. A structured decision-making approach that addresses regulatory requirements, operational needs, and sustainability goals empowers pharmaceutical manufacturers to make informed, responsible glove selections, optimizing performance and supporting broader corporate sustainability initiatives.

Walid El Azab is a freelance writer commissioned by O&M Halyard to write this article and operates independently of the company

Reference:

1. REGULATION (EU) 2016/425 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 9 March 2016 on personal protective equipment and repealing Council Directive 89/686/EEC : <https://eur-lex.europa.eu/legal-content/EN/TXT/HTML/?uri=CELEX:32016R0425>
2. Rego A, Roley L. In-use barrier integrity of gloves: latex and nitrile superior to vinyl. *Am J Infect Control*. 1999 Oct;27(5):405-10. doi: 10.1016/s0196-6553(99)70006-4. Erratum in: *Am J Infect Control* 1999 Dec;27(6):542. PMID: 10511487.
3. Phalen RN, Le T, Wong WK. Changes in chemical permeation of disposable latex, nitrile, and vinyl gloves exposed to simulated movement. *J Occup Environ Hyg*. 2014;11(11):716-21. doi: 10.1080/15459624.2014.908259. PMID: 24689368; PMCID: PMC4177299.
4. JO SAWYER, ALLAN BENNETT, Comparing the Level of Dexterity offered by Latex and Nitrile SafeSkin Gloves, *The Annals of Occupational Hygiene*, Volume 50, Issue 3, April 2006, Pages 289–296, <https://doi.org/10.1093/annhyg/mei066>
5. Chadwick, R. Comparison of nitrile and latex examination gloves. *Br Dent J* 196, 685 (2004). <https://doi.org/10.1038/sj.bdj.4811350>
6. Slater JE. Latex allergy. *J Allergy Clin Immunol*. 1994 Aug;94(2 Pt 1):139-49; quiz 150. doi: 10.1016/0091-6749(94)90031-0. PMID: 8064066.
7. Hunt LW, Fransway AF, Reed CE, Miller LK, Jones RT, Swanson MC, Yunginger JW. An epidemic of occupational allergy to latex involving health care workers. *J Occup Environ Med*. 1995 Oct;37(10):1204-9. doi: 10.1097/00043764-199510000-00011. PMID: 8542340.
8. Heese A, van Hintzenstern J, Peters KP, Koch HU, Hornstein OP. Allergic and irritant reactions to rubber gloves in medical health services. Spectrum, diagnostic approach, and therapy. *J Am Acad Dermatol*. 1991 Nov;25(5 Pt 1):831-9. doi: 10.1016/s0190-9622(08)80977-2. Erratum in: *J Am Acad Dermatol* 1992 Mar;26(3 Pt 2):403. PMID: 1802906.
9. Hansen A, Brans R, Sonnsman F. Allergic contact dermatitis to rubber accelerators in protective gloves: Problems, challenges, and solutions for occupational skin protection. *Allergol Select*. 2021 Oct 22;5:335-344. doi: 10.5414/ALX02265E. PMID: 34734158; PMCID: PMC8561832.
10. EudraLex The Rules Governing Medicinal Products in the European Union, Volume 4: EU Guidelines to Good Manufacturing Practice Medicinal Products for Human and Veterinary Use, Annex 1 Manufacture of Sterile Medicinal Products, (2024).
11. Esmizadeh, E., et al. (2021). "Stability of nitrile and vinyl latex gloves under repeated disinfection cycles." **Materials Today Sustainability**, 11-12: 100067.
12. Phalen, R. N., et al. (2014). "Changes in Chemical Permeation of Disposable Latex, Nitrile, and Vinyl Gloves Exposed to Simulated Movement." **Journal of Occupational and Environmental Hygiene**, 11(11): 716-721.
13. Gao P, Horvatin M, Niezgoda G, Weible R, Shaffer R. Effect of multiple alcohol-based hand rub applications on the tensile properties of thirteen brands of medical exam nitrile and latex gloves. *J Occup Environ Hyg*. 2016 Dec;13(12):905-914. doi: 10.1080/15459624.2016.1191640. PMID: 27224677; PMCID: PMC7157957.
14. Sandle T., Glove disinfection and aseptic technique: Creating a schema for the cleanroom and laboratory, *EJPPS EUROPEAN JOURNAL OF PARENTERAL AND PHARMACEUTICAL SCIENCES* · July 2023
15. Amgen to Achieve Carbon Neutrality by 2027, accessed on Dec 24: <https://www.amgen.com/newsroom/press-releases/2021/01/amgen-to-achieve-carbon-neutrality-by-2027#:~:text=Amgen%20will%20invest%20more%20than,efficiencics%20over%20the%20same%20period> and <https://www.amgen.com/responsibility/healthy-planet/environmental-sustainability/approach>
16. Novartis Environmental Sustainability Strategy, accessed on Dec 2024: https://www.novartis.com/sites/novartis_com/files/novartis-environmental-sustainability-strategy-objectives.pdf
17. Novo Nordisk, Zero Environmental impact; accessed on Dec 2024: https://www.novonordisk.com/sustainable-business/zero-environmental-impact.html?utm_source=chatgpt.com
18. Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste and repealing certain Directives
19. United States Environmental Protection Agency, Law & Executive order: 42 U.S.C. §6901 et seq. (1976) Summary of the Resource Conservation and Recovery Act visited on Dec 2024: <https://www.epa.gov/laws-regulations/summary-resource-conservation-and-recovery-act>
20. European Environment Agency, measuring the quality of recycling, accessed on Dec 2024: <https://www.eea.europa.eu/en/analysis/publications/measuring-the-quality-of-recycling>
21. The Climate Crisis: Working together for Future Generations, accessed on Dec 2024: <https://www.state.gov/policy-issues/climate-crisis/>
22. Directive (EU) 2022/2464 of the European Parliament and of the Council of 14 December 2022 amending Regulation (EU) No 537/2014, Directive 2004/109/EC, Directive 2006/43/EC and Directive 2013/34/EU, as regards corporate sustainability reporting
23. Deep Dive: The EU Corporate Sustainability Due Diligence Directive & Implications for Companies; Accessed on Jan 2025: <https://fiscalnote.com/blog/eu-corporate-sustainability-due-diligence-directive>

24. US Department of Labor, Wage and Hour division, Wages and the Fair Labor Standard Act, accessed on Jan 205: <https://www.dol.gov/agencies/whd/flsa>
25. Responsible Glove Alliance, accessed on Jan 2025: <https://www.responsibleglove.org/about/rga-overview/>
26. Hansen A, Brans R, Sonnsmann F. Allergic contact dermatitis to rubber accelerators in protective gloves: Problems, challenges, and solutions for occupational skin protection. Allergol Select. 2021 Oct 22;5:335-344. doi: 10.5414/ALX02265E. PMID: 34734158; PMCID: PMC8561832.
27. Directive (EU) 2024/1760 of the European Parliament and of the Council of 13 June 2024 on corporate sustainability due diligence and amending Directive (EU) 2019/1937 and Regulation (EU) 2023/2859
28. SEC Adopts Rules to Enhance and Standardize Climate-Related Disclosures for Investors, Accessed on Jan 2025 <https://www.sec.gov/newsroom/press-releases/2024-31>

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Walid El Azab

QP Pro Services & QPM Consulting
Co-Fondateur et Consultant Senior

Trusted partner for Compliance,
Continuity and Business Excellence

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Walid El Azab – QP Pro Services & QPM Consulting
Co-Fondateur et Consultant Senior

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Référence:

1. REGULATION (EU) 2016/425 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 9 March 2016 on personal protective equipment and repealing Council Directive 89/686/EEC : <https://eur-lex.europa.eu/legal-content/EN/TXT/HTML/?uri=CELEX:32016R0425>
2. Rego A, Roley L. In-use barrier integrity of gloves: latex and nitrile superior to vinyl. *Am J Infect Control*. 1999 Oct;27(5):405-10. doi: 10.1016/s0196-6553(99)70006-4. Erratum in: *Am J Infect Control* 1999 Dec;27(6):542. PMID: 10511487.
3. Phalen RN, Le T, Wong WK. Changes in chemical permeation of disposable latex, nitrile, and vinyl gloves exposed to simulated movement. *J Occup Environ Hyg*. 2014;11(11):716-21. doi: 10.1080/15459624.2014.908259. PMID: 24689368; PMCID: PMC4177299.
4. JO SAWYER, ALLAN BENNETT, Comparing the Level of Dexterity offered by Latex and Nitrile SafeSkin Gloves, *The Annals of Occupational Hygiene*, Volume 50, Issue 3, April 2006, Pages 289–296, <https://doi.org/10.1093/annhyg/mei066>
5. Chadwick, R. Comparison of nitrile and latex examination gloves. *Br Dent J* 196, 685 (2004). <https://doi.org/10.1038/sj.bdj.4811350>
6. Slater JE. Latex allergy. *J Allergy Clin Immunol*. 1994 Aug;94(2 Pt 1):139-49; quiz 150. doi: 10.1016/0091-6749(94)90031-0. PMID: 8064066.
7. Hunt LW, Fransway AF, Reed CE, Miller LK, Jones RT, Swanson MC, Yunginger JW. An epidemic of occupational allergy to latex involving health care workers. *J Occup Environ Med*. 1995 Oct;37(10):1204-9. doi: 10.1097/00043764-199510000-00011. PMID: 8542340.
8. Heese A, van Hintzenstern J, Peters KP, Koch HU, Hornstein OP. Allergic and irritant reactions to rubber gloves in medical health services. Spectrum, diagnostic approach, and therapy. *J Am Acad Dermatol*. 1991 Nov;25(5 Pt 1):831-9. doi: 10.1016/s0190-9622(08)80977-2. Erratum in: *J Am Acad Dermatol* 1992 Mar;26(3 Pt 2):403. PMID: 1802906.
9. Hansen A, Brans R, Sonsmann F. Allergic contact dermatitis to rubber accelerators in protective gloves: Problems, challenges, and solutions for occupational skin protection. *Allergol Select*. 2021 Oct 22;5:335-344. doi: 10.5414/ALX02265E. PMID: 34734158; PMCID: PMC8561832.
10. EudraLex The Rules Governing Medicinal Products in the European Union, Volume 4: EU Guidelines to Good Manufacturing Practice Medicinal Products for Human and Veterinary Use, Annex 1 Manufacture of Sterile Medicinal Products, (2024).
11. Esmizadeh, E., et al. (2021). "Stability of nitrile and vinyl latex gloves under repeated disinfection cycles." *Materials Today Sustainability*, 11-12: 100067.
12. Phalen, R. N., et al. (2014). "Changes in Chemical Permeation of Disposable Latex, Nitrile, and Vinyl Gloves Exposed to Simulated Movement." *Journal of Occupational and Environmental Hygiene*, 11(11): 716-721.
13. Gao P, Horvatin M, Niezgoda G, Weible R, Shaffer R. Effect of multiple alcohol-based hand rub applications on the tensile properties of thirteen brands of medical exam nitrile and latex gloves. *J Occup Environ Hyg*. 2016 Dec;13(12):905-914. doi: 10.1080/15459624.2016.1191640. PMID: 27224677; PMCID: PMC7157957.
14. Sandle T., Glove disinfection and aseptic technique: Creating a schema for the cleanroom and laboratory, *EJPPS EUROPEAN JOURNAL OF PARENTERAL AND PHARMACEUTICAL SCIENCES* · July 2023
15. Amgen to Achieve Carbon Neutrality by 2027, accessed on Dec 24: <https://www.amgen.com/newsroom/press-releases/2021/01/amgen-to-achieve-carbon-neutrality-by-2027#:~:text=Amgen%20will%20invest%20more%20than,efficiencies%20over%20the%20same%20period> and <https://www.amgen.com/responsibility/healthy-planet/environmental-sustainability/approach>
16. Novartis Environmental Sustainability Strategy, accessed on Dec 2024: <https://www.novartis.com/sites/novartis.com/files/novartis-environmental-sustainability-strategy-objectives.pdf>
17. Novo Nordisk, Zero Environmental impact; accessed on Dec 2024: https://www.novonordisk.com/sustainable-business/zero-environmental-impact.html?utm_source=chatgpt.com
18. Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste and repealing certain Directives
19. United States Environmental Protection Agency, Law & Executive order: 42 U.S.C. §6901 et seq. (1976) Summary of the Resource Conservation and Recovery Act visited on Dec 2024: <https://www.epa.gov/laws-regulations/summary-resource-conservation-and-recovery-act>
20. European Environment Agency, measuring the quality of recycling, accessed on Dec 2024: <https://www.eea.europa.eu/en/analysis/publications/measuring-the-quality-of-recycling>
21. The Climate Crisis: Working together for Future Generations, accessed on Dec 2024: <https://www.state.gov/policy-issues/climate-crisis/>
22. Directive (EU) 2022/2464 of the European Parliament and of the Council of 14 December 2022 amending Regulation (EU) No 537/2014, Directive 2004/109/EC, Directive 2006/43/EC and Directive 2013/34/EU, as regards corporate sustainability reporting
23. Deep Dive: The EU Corporate Sustainability Due Diligence Directive & Implications for Companies; Accessed on Jan 2025: <https://fiscalnote.com/blog/eu-corporate-sustainability-due-diligence-directive>
24. US Department of Labor, Wage and Hour division, Wages and the Fair Labor Standard Act, accessed on Jan 205: <https://www.dol.gov/agencies/whd/flsa>
25. Responsible Glove Alliance, accessed on Jan 2025: <https://www.responsibleglove.org/about/rga-overview/>
26. Hansen A, Brans R, Sonsmann F. Allergic contact dermatitis to rubber accelerators in protective gloves: Problems, challenges, and solutions for occupational skin protection. *Allergol Select*. 2021 Oct 22;5:335-344. doi: 10.5414/ALX02265E. PMID: 34734158; PMCID: PMC8561832.
27. Directive (EU) 2024/1760 of the European Parliament and of the Council of 13 June 2024 on corporate sustainability due diligence and amending Directive (EU) 2019/1937 and Regulation (EU) 2023/2859
28. SEC Adopts Rules to Enhance and Standardize Climate-Related Disclosures for Investors, Accessed on Jan 2025 <https://www.sec.gov/newsroom/press-releases/2024-31>

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SHAPING COMPLIANCE & EXCELLENCE WITH **HIGHLY QUALIFIED** & **PRAGMATIC EXPERTS**

Who We Are

Who We Are QP Pro Services is an international, highly specialized consultancy firm with over 20 years of experience. Led by a team of Highly Qualified Persons (QPs) and Responsible Persons (RPs), we are dedicated to supporting the pharmaceutical and life sciences industries.

WHY CHOOSE US?

- Achieving operational excellence and regulatory compliance
- Empowering Decision
- Making with Data
- Driving Innovation & Improvement Growth
- Fostering Cost-Effective Excellence
- Securing Long-lasting Partnerships
- Providing Next Gen Experts: Highly skilled QP network to meet the need or lack of the next generation of QP's

Our Mission

we strive to be the trusted partner in compliance and operational excellence. Our mission is to deliver pragmatic and tailored solutions that empower the pharmaceutical and life sciences industries to achieve regulatory compliance, drive innovation, and ensure product integrity at every stage.

OUR CORE COMPETENCES

Customer Focused
Teamwork / Collaboration
Deliver & Execute results / Striving For excellence
Pragmatism Leadership (self & people)
Continuous Development
Communication & Change Management
Emotional Intelligence
Decision Making
Job Knowledge/Business Acumen

Our values

EXCELLENCE
COLLABORATION
INTEGRITY
CUSTOMER FIRST
ACCOUNTABILITY

Contact us Now

info@qpproservices.com

www.qpproservices.com

Our Social Media

[linkedin|qpproservices](#)

THE REFERENCE IN PHARMACEUTICAL COMPLIANCE AND OPERATIONAL EXCELLENCE

We ensure tailored services for
the life sciences industry.

CONTACT US

Stay Ahead with QP Pro Services

Take full advantage of our events, publications, and training opportunities to stay connected and enhance your expertise. By leveraging data and digital tools, we empower QPs and RPs to make pragmatic decisions that balance compliance with business growth. Together, we ensure your organization remains a leader in pharmaceutical excellence.

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SERVICES

DRUG PRODUCT EU IMPORTATION

End-to-End EU importation,
distribution and product
release.

A network of high skilled QP's,
RP's and Partners.

CONSULTING & OPERATIONAL SUPPORT

Tailored training and
consultancy in quality
assurance, contamination
control, sterility assurance,
manufacturing, validation,
auditing, inspection readiness,
regulatory affairs,
pharmacovigilance, and data-
driven solutions or ad-interim
management.

QP & RP

Batch certification, distribution
for clinical and commercial
phases, including post-marketing
and pharmacovigilance.

Ensuring QP/RP back-up or
contingency solutions to
safeguard product integrity and
compliance at every stage.